

Posterior Pole Asymmetry Analysis



Asymmetry is a hallmark of glaucoma. Posterior Pole Asymmetry Analysis can help identify early glaucomatous damage.

Dr. Sanjay Asrani
(Duke University)



Glaucoma is a disease that manifests by loss of ganglion cells and axons across the central posterior pole, where the ganglion cells are most concentrated. Retinal thickness reflects glaucomatous damage by marked thinning in the zone surrounding the fovea, extending towards the optic nerve head.

The SPECTRALIS OCT Posterior Pole Asymmetry Analysis combines mapping of the posterior pole retinal thickness with asymmetry analysis between eyes and between hemispheres of each eye.

RNFL measurements combined with retina thickness measurement gives a much more robust parameter for glaucoma.



Case 1:
Retinal thickness map of a glaucomatous eye and its fellow eye. Note the severe localized thinning in the inferotemporal region of the right eye. Asymmetry between eyes and hemispheres illustrates damage.



Case 2:
Retinal thickness map of the right eye shows a significant thinning infero- and superotemporally. A clear asymmetry between eyes becomes visible.

How to Interpret the Asymmetry Analysis

Posterior Pole Asymmetry Single Exam Report OU

SPECTRALIS® Tracking Laser Tomography

Patient: Posterior Pole, - Case 4
 Patient ID: ...
 Diagnosis: ...

DOB: 29/Aug/1949
 Exam: 30/Nov/2009
 Comment: ...

Gender: F



Posterior Pole Retina Thickness Map
 Displays the retinal thickness over the entire posterior pole (30° x 25° OCT volume scan) for each eye.





Compressed Color Scale
 Used to localize even the smallest differences in retinal thickness.

8x8 Analysis Grid
 An 8x8 grid is positioned symmetrically to the fovea-disc axis. For each cell of the grid the mean retinal thickness is given.

Hemisphere Analysis
 Displays the asymmetry between the superior and inferior hemisphere. The fovea-disc axis is the horizontal symmetry line. For each cell of one hemisphere, the mean retina thickness is compared to the value in the corresponding cell for the opposite hemisphere.

I-S Asymmetry Map
 The lower half compares the inferior to the superior hemisphere.

OD



OS

OU Asymmetry
 The mean retinal thickness in each cell of the 8x8 grid in one eye is compared to the thickness in the corresponding cell of the fellow eye.

OS-OD Asymmetry Map
 Compares the left eye to the right eye.

Asymmetry Color Scale
 Darker grey indicating larger differences. The closer the value is to zero (white color), the better the symmetry.

Mean Thickness
 Represents the mean retinal thickness for the superior and inferior hemisphere as well as the total mean thickness over the entire 8x8 grid.

Hemisphere Asymmetry S-I		Mean Thickness		Hemisphere Asymmetry S-I	
	Superior (S) 299 µm Inferior (I) 272 µm	Superior (S) 299 µm Inferior (I) 295 µm			
Total 285 µm		Total 297 µm			

Notes:
 Date: 20/08/2010
 Signature:
 Software Version: 5.3.0.22
 www.HeidelbergEngineering.com
 Posterior Pole Asymmetry Single Exam Report OU



Headquarters
 Heidelberg Engineering GmbH
 Tiegartenstr. 15
 69121 Heidelberg · Germany
 Tel: +49 6221 6463-0 · Fax +49 6221 6463-62
 www.HeidelbergEngineering.com

USA
 Heidelberg Engineering, Inc.
 1499 Poinsettia Avenue, Suite 160
 Vista, CA 92081
 Tel: +1 760-598-3770 · Fax +1 760-598-3060
 www.HeidelbergEngineering.com