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1. KIM KY et al.:Anterior segment characteristics in normal and keratoconus eyes evaluated with a new type of swept-source optical coherence tomography.PLoS One 2022; 17(9):e0274071.DOI:10.1371/journal.pone.0274071.
2. HERBER R et al.:Agreement and repeatability of corneal tomography in healthy eyes using a new swept-source OCT a rotating Scheimpflug camera and a dual Scheimpflug-Placido system.J Cataract Refract Surg 2022; 48(2):190-198.
3. KIM K et al.:Diagnostic Validation of the Screening Corneal Objective Risk of Ectasia Analyzer Evaluated by Swept Source Optical Coherence Tomography for Keratoconus in an Asian Population.Bioengineering (Basel) 2023; 10(11):1335.
4. The Heidelberg Engineering ANTERION Anterior Segment Cornea and IOL Precision and Agreement Study [data on file].
5. TAÑÁ-RIVERO P et al.:Repeatability of whole-cornea measurements using a new swept-source optical coherence tomographer.Eur J Ophthalmol 2021; 31(4):1709-1719.
6. ESCOLANO SERRANO J et al.:Intraobserver Repeatability of Tomographic Pachymetric and Anatomical Measurements in Healthy Eyes Using a New Swept-Source Optical Coherence Topographer.Cornea 2022; 41(5):598-603.
7. HERBER R et al.:Comparison of corneal tomography using a novel swept-source optical coherence tomographer and rotating Scheimpflug system in normal and keratoconus eyes: repeatability and agreement analysis.Eye Vis (Lond) 2022; 9(1):19.
8. PÉREZ-BARTOLOMÉ F et al.:Anterior-Segment Swept-Source Ocular Coherence Tomography and Scheimpflug Imaging Agreement for Keratometry and Pupil Measurements in Healthy Eyes.J Clin Med 2021; 10(24):5789.
9. SAAD A et al.:Discrimination between keratoconus forme fruste keratoconus and normal eyes using a novel OCT-based tomographer.J Cataract Refract Surg 2023; 49(11):1092-1097.
10. SCHIANO-LOMORIELLO D et al.:Repeatability of automated measurements by a new anterior segment optical coherence tomographer and biometer and agreement with standard devices.Sci Rep 2021; 11(1):983.
11. CHENG SM et al.:Repeatability of a new swept-source optical coherence tomographer and agreement with other three optical biometers.Graefes Arch Clin Exp Ophthalmol 2022; 260(7):2271-2281.
12. TAÑÁ-RIVERO P et al.:Agreement between 2 swept-source OCT biometers and a Scheimpflug partial coherence interferometer.J Cataract Refract Surg 2021; 47(4):488-495.
13. TAÑÁ-RIVERO P et al.:Lens-vault analysis and its correlation with other biometric parameters using swept-source OCT.J Optom 2022; 15(1):88-99.
14. TANA-SANZ P et al.:Agreement of predicted intraocular lens power using swept-source optical coherence tomography and partial coherence interferometry.Expert Rev Med Devices 2021; 18(12):1219-1234.
15. RUIZ-MESA R et al.:Ocular biometric repeatability using a new high-resolution swept-source optical coherence tomographer.Expert Rev Med Devices 2020; 17(6):591-597.
16. CHENG SM et al.:Repeatability and Agreement of Two Swept-Source Optical Coherence Tomographers for Anterior Segment Parameter Measurements.J Glaucoma 2022; 31(7):602-608.
17. XIE X et al.:Age- and refraction-related changes in anterior segment anatomical structures measured by swept-source anterior segment OCT.PLoS One 2020; 15(10):e0240110. DOI:10.1371/journal.pone.0240110.
18. PARDESHI AA et al.:Intradevice Repeatability and Interdevice Agreement of Ocular Biometric Measurements:A Comparison of Two Swept-Source Anterior Segment OCT Devices.Transl Vis Sci Technol 2020; 9(9):14.
19. CHAN PP et al.:Anterior chamber angle imaging with swept-source optical coherence tomography:Comparison between CASIAII and ANTERION.Sci Rep 2020; 10(1):18771.
20. SAVINI G:Posterior Corneal Measurements by 4 Devices.IOL Power Club; November 2023; Palm Springs (not published).
21. PÉREZ-BARTOLOMÉ F et al.:Agreement between anterior segment swept source-OCT and Scheimpflug imaging corneal aberration measurements in healthy eyes.Eur J Ophthalmol 2022; 32(6):3363-3371.
22. YUN JS et al.:Evaluation of angle-to-angle and spur-to-spur using swept source optical coherence tomography in different refractive error.PLoS One 2022; 17(11):e0277703.
23. MONTÉS-MICÓ R et al.:Angle-to-angle and spur-to-spur distance analysis with high-resolution optical coherence tomography.Eye Vis (Lond) 2020; 7:42.
24. MONTÉS-MICÓ R et al.:Assessment of anterior segment measurements using a high-resolution imaging device. Expert Rev Med Devices 2020; 17(9):969-979.
25. TAÑÁ-RIVERO P et al.:Agreement of white-to-white measurements with swept-source OCT Scheimpflug and color LED devices.Int Ophthalmol 2021; 41(1):57-65.
26. LIN F et al.:Comparison of visual outcomes and higher-order aberrations between FS-LASIK and SMI-LIKE for moderate to high hyperopia: a 2-year result.Cornea 2022; 10:1097.
27. DE OCA IM et al.:IOL calculations in short, long, and postrefractive eyes.International Ophthalmology Clinics 2016; 56(3):49-70.
28. HAIGIS W:Preoperative intraocular lens power calculations in problem eyes.Zeitschrift fur Medizinische Physik 2007; 17(1):45-54.
29. OLSEN T:Calculation of intraocular lens power: a review.Acta Ophthalmologica Scandinavica 2007; 85(5):472-485.
30. KOTHARI SS et al.:Recent developments in the intraocular lens formulae:An update.Seminars in Ophthalmology 2023; 38(2):143-150.
31. CHEN Y et al.:Silicone oil in vitreoretinal surgery: indications, complications, new developments and alternative long-term tamponade agents.Acta ophthalmologica 2021; 99(3):240-250.
32. KANCLERZ P, GRZYBOWSKI A:Complications associated with the use of expandable gases in vitrectomy.Journal of ophthalmology 2018:8606494.

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Heidelberg Engineering GmbH, 德国海德堡公司 · Max-Jarecki-Str.8 · 69115 Heidelberg · Germany (德国)  
电话号码: +49 6221 64630 · 传真号: +49 6221 646362

**AUS**

Heidelberg Engineering Pty Ltd · Suite E5, 63 Turner Street · Port Melbourne 3207 · Victoria  
电话号码: +61 396 392 125 · 传真号: +61 396 392 127

**CH**

Heidelberg Engineering GmbH, 德国海德堡公司 · Schulstrasse 161 · 8105 Regensdorf  
电话号码: +41 44 8887 020 · 传真号: +41 44 8887 024

**FIN**

Heidelberg Engineering GmbH, 德国海德堡公司 · Suomen sivuliike · Fönsenintie 1 · 00370 Helsinki  
电话号码: +358 505 226 963

**UK**

Heidelberg Engineering Ltd. · 55 Marlowes · Hemel Hempstead · Hertfordshire HP1 1LE  
电话号码: +44 1442 502 330 · 传真号: +44 1442 242 386

<http://www.HeidelbergEngineering.com>