



JANE, DOE

ID: 0005446181
Sex: F
Birthdate: 16.03.1949

Exam Date: 17.01.2022
Exam Description: External Patient:



Advanced OCT Report



JANE, DOE

ID: 0005446181
Sex: F
Birthdate: 16.03.1949

Exam Date: 17.01.2022
Exam Description: External Patient:

This report cannot be considered a clinical diagnosis report. AI/ML components/models/applications used for the report are intended for research purposes and decision support assistance only.




Scan severity


Right (OD)

Severity: Red

64 / 128



No Severity (green): 0 / 128
Moderate Severity (yellow): 17 / 128
Severe (red): 111 / 128




Retina layers


Right (OD)

64 / 128


- RNFL (Retinal nerve fiber layer) ✓
- GCL (Ganglion cell layer) ✓
- IPL (Inner plexiform layer) ✓
- INL (Inner nuclear layer) ✓
- OPL (Outer plexiform layer) ✓
- ONL (Outer nuclear layer + Henle`s nerve fiber layer) ✓
- ELM (External limiting membrane) ✓
- MZ (Myoid zone of the photoreceptors) ✓
- EZ (Ellipsoid zone of the photoreceptors) ✓
- OS (Outer segments of the photoreceptors) ✓
- RPE (Retinal pigment epithelium) ✓
- BM (Bruch`s membrane) ✓




ETDRS Calculations



1/3/6 mm average thickness, μm



GCC thickness (RNFL+GCL+IPL)



GCC thickness (GCL+IPL)

JANE, DOE

ID: 0005446181
Sex: F
Birthdate: 16.03.1949

Exam Date: 17.01.2022
Exam Description: External Patient:

This report cannot be considered a clinical diagnosis report. AI/ML components/models/applications used for the report are intended for research purposes and decision support assistance only.




Pathologies detection

Right (OD)

Segmentation

64 / 128

- Diffuse edema ✓
- Disorganization of retinal inner layers ✗
- Epiretinal fibrosis ✓
- Intraretinal cystoid fluid ✓
- Intraretinal hyperreflective foci ✓
- Subretinal fluid ✓



Classification:

Diabetic macular edema


Diabetic retinopathy

Epiretinal fibrosis

Segmentation

23 / 128

- Diffuse edema ✓
- Epiretinal fibrosis ✓
- Hard exudates ✓
- Hyperreflectivity of neurosensory retina ✓
- Intraretinal hyperreflective foci ✓
- RPE disruption ✓
- Shadowing ✓



Classification:

Diabetic macular edema

Diabetic retinopathy

JANE, DOE

ID: 0005446181
Sex: F
Birthdate: 16.03.1949

Exam Date: 17.01.2022
Exam Description: External Patient:

This report cannot be considered a clinical diagnosis report. AI/ML components/models/applications used for the report are intended for research purposes and decision support assistance only.




Pathologies detection

Right (OD)

Segmentation

17 / 128

- Diffuse edema ✓
- Ellipsoid zone disruption ✓
- Epiretinal fibrosis ✓
- Hard exudates ✓
- Hyperreflectivity of neurosensory retina ✓
- Intraretinal hyperreflective foci ✓
- RPE disruption ✓
- Retinal wrinkling ✓
- Shadowing ✓




Classification: Diabetic macular edema Diabetic retinopathy

Segmentation

25 / 128

- Diffuse edema ✓
- Ellipsoid zone disruption ✓
- Epiretinal fibrosis ✓
- Hard exudates ✓
- Hyperreflectivity of neurosensory retina ✓
- Intraretinal hyperreflective foci ✓
- RPE disruption ✓
- Shadowing ✓



Classification: Diabetic macular edema Diabetic retinopathy

JANE, DOE

ID: 0005446181

Sex: F

Birthdate: 16.03.1949

Exam Date: 17.01.2022

Exam Description: External Patient:

This report cannot be considered a clinical diagnosis report. AI/ML components/models/applications used for the report are intended for research purposes and decision support assistance only.



Pathologies Area and Volume

Right (OD)

Scan: 64

Name	Volume (mm ³) (ETDRS 1 mm)	Volume (mm ³) (ETDRS 3 mm)	Volume (mm ³) (ETDRS 6 mm)	B-scan area (mm ²)	B-scan area view	En face area view (mm ² * 1000)
Diffuse edema	0.3026	2.8888	10.5276	2.4497		
Disorganization of retinal inner layers	0.0032	0.0521	0.1334	0.1868		
Epiretinal fibrosis	0.0018	0.0591	0.2004	0.0322		
Intraretinal cystoid fluid	0.0000	0.0051	0.0188	0.0047		
Intraretinal hyperreflective foci	0.0006	0.0058	0.0375	0.0181		
Subretinal fluid	0.0000	0.0057	0.0057	0.0076		

JANE, DOE

ID: 0005446181
 Sex: F
 Birthdate: 16.03.1949

Exam Date: 17.01.2022
 Exam Description: External Patient:

This report cannot be considered a clinical diagnosis report. AI/ML components/models/applications used for the report are intended for research purposes and decision support assistance only.



Pathologies Area and Volume Right (OD)

Scan: 23

Name	Volume (mm ³) (ETDRS 1 mm)	Volume (mm ³) (ETDRS 3 mm)	Volume (mm ³) (ETDRS 6 mm)	B-scan area (mm ²)	B-scan area view	En face area view (mm ² * 1000)
Diffuse edema	0.3026	2.8888	10.5276	1.9558		
Epiretinal fibrosis	0.0018	0.0591	0.2004	0.0495		
Hard exudates	0.0000	0.0003	0.0394	0.0539		
Hyperreflectivity of neurosensory retina	0.0000	0.0197	0.2015	0.1787		
Intraretinal hyperreflective foci	0.0006	0.0058	0.0375	0.0206		
RPE disruption	0.0053	0.0141	0.0265	0.0291		
Shadowing	0.0000	0.0000	0.1250	0.1263		

JANE, DOE

ID: 0005446181
 Sex: F
 Birthdate: 16.03.1949

Exam Date: 17.01.2022
 Exam Description: External Patient:

This report cannot be considered a clinical diagnosis report. AI/ML components/models/applications used for the report are intended for research purposes and decision support assistance only.



Pathologies Area and Volume Right (OD)

Scan: 17

Name	Volume (mm ³) (ETDRS 1 mm)	Volume (mm ³) (ETDRS 3 mm)	Volume (mm ³) (ETDRS 6 mm)	B-scan area (mm ²)	B-scan area view	En face area view (mm ² * 1000)
Diffuse edema	0.3026	2.8888	10.5276	1.825		
Epiretinal fibrosis	0.0018	0.0591	0.2004	0.0222		
Hard exudates	0.0000	0.0003	0.0394	0.0092		
Hyperreflectivity of neurosensory retina	0.0000	0.0197	0.2015	0.2677		
Intraretinal hyperreflective foci	0.0006	0.0058	0.0375	0.0311		
Ellipsoid zone disruption	0.0000	0.0008	0.0011	0.0044		
Retinal wrinkling	0.0000	0.0006	0.0020	0.0033		

JANE, DOE

ID: 0005446181
Sex: F
Birthdate: 16.03.1949



Exam Date: 17.01.2022
Exam Description: External Patient:

This report cannot be considered a clinical diagnosis report. AI/ML components/models/applications used for the report are intended for research purposes and decision support assistance only.



Pathologies Area and Volume Right (OD)

Scan: 17

Name	Volume (mm ³) (ETDRS 1 mm)	Volume (mm ³) (ETDRS 3 mm)	Volume (mm ³) (ETDRS 6 mm)	B-scan area (mm ²)	B-scan area view	En face area view (mm ² * 1000)
● RPE disruption	0.0053	0.0141	0.0265	0.0262		
● Shadowing	0.0000	0.0000	0.1250	0.0343		

JANE, DOE

ID: 0005446181

Sex: F

Birthdate: 16.03.1949

Exam Date: 17.01.2022

Exam Description: External Patient:

This report cannot be considered a clinical diagnosis report. AI/ML components/models/applications used for the report are intended for research purposes and decision support assistance only.



Pathologies Area and Volume

Right (OD)

Scan: 25

Name	Volume (mm ³) (ETDRS 1 mm)	Volume (mm ³) (ETDRS 3 mm)	Volume (mm ³) (ETDRS 6 mm)	B-scan area (mm ²)	B-scan area view	En face area view (mm ² * 1000)
Diffuse edema	0.3026	2.8888	10.5276	1.8671		
Epi-retinal fibrosis	0.0018	0.0591	0.2004	0.0502		
Hard exudates	0.0000	0.0003	0.0394	0.0484		
Hyperreflectivity of neurosensory retina	0.0000	0.0197	0.2015	0.035		
Intraretinal hyperreflective foci	0.0006	0.0058	0.0375	0.0165		
Ellipsoid zone disruption	0.0000	0.0008	0.0011	0.0061		
RPE disruption	0.0053	0.0141	0.0265	0.0308		

JANE, DOE

ID: 0005446181
Sex: F
Birthdate: 16.03.1949

Exam Date: 17.01.2022
Exam Description: External Patient:

This report cannot be considered a clinical diagnosis report. AI/ML components/models/applications used for the report are intended for research purposes and decision support assistance only.



Pathologies Area and Volume Right (OD)

Scan: 25

Name	Volume (mm ³) (ETDRS 1 mm)	Volume (mm ³) (ETDRS 3 mm)	Volume (mm ³) (ETDRS 6 mm)	B-scan area (mm ²)	B-scan area view	En face area view (mm ² * 1000)
● Shadowing	0.0000	0.0000	0.1250	0.1228		

JANE, DOE

ID: 0005446181
Sex: F
Birthdate: 16.03.1949

Exam Date: 17.01.2022
Exam Description: External Patient:

This report cannot be considered a clinical diagnosis report. AI/ML components/models/applications used for the report are intended for research purposes and decision support assistance only.



Glaucoma Risk analysis is unavailable for this examination.
This could be due to:

Pathologies present in the scans affecting the GCC complex, making accurate calculation impossible.
One of the thicknesses out of boundaries of 125 μm