



# Fluid Monitor Report

nAMD (neovascular age-related macular degeneration)



Heidelberg Engineering / Spectralis / OS	Previous Report	Current Report
Visit	2022-08-13 00:00 (2 months and 3 days prior)	2022-10-16 00:00
Scan Pattern	6.0mm x 6.0mm, 97 B-scans, distance: 60µm	6.0mm x 6.0mm, 97 B-scans, distance: 60µm

Compartment		Measured Volume (in Nanoliter)			
		ETDRS 1mm		ETDRS 6mm	
		Previous Report	Current Report	Previous Report	Current Report
<b>IRF</b>	intraretinal fluid	10.00	3.00 ↓ 7.00	10.00	3.00 ↓ 7.00
<b>SRF</b>	subretinal fluid	18.00	5.00 ↓ 13.00	18.00	5.00 ↓ 13.00
<b>PED</b>	pigment epithelial detachment	23.00	7.00 ↓ 16.00	23.00	7.00 ↓ 16.00

## Progression in 1mm Region

